













1. Prepare injection site per hospital protocol. Ensure skin is clean and dry.
2. Remove HaloGUARD™ Protective Disc with CHG from package using sterile technique.
3. Place HaloGUARD™ Protective Disc with CHG around compatible device. Ensure print side faces upward.
4. Ensure compatible device rests upon slight portion of HaloGUARD™ Protective Disc with CHG. Do not compress disc with compatible device connector. The disc must remain in contact with skin to maximize efficacy.
5. Secure compatible device and HaloGUARD™ Protective Disc with CHG to skin, using a bioadhesive dressing. Ensure complete skin contact.
6. Change HaloGUARD™ Protective Disc with CHG according to hospital protocol.
7. Antibiogram of HaloGUARD™ Protective Disc with CHG was evaluated for up to 7 days in-vitro, therefore, dressing changes should occur at least once every 7 days. Highly exuding wounds require more frequent dressing changes.

**INSTRUCTIONS FOR USE**  
(Read Carefully Before Use)

|   |                                     |
|---|-------------------------------------|
| Labelling Symbols   |                                     |
|  | Do not re-use                       |
|  | Caution! See instructions for use   |
|  | Do not use if package is damaged    |
|  | Do not re-sterilize using radiation |
|  | Temperature limits                  |
|  | STERILE R                           |
|  | QTY Quantity                        |
|  | Use by Date                         |
|  | Manufacturer                        |
|  | REF Catalogue number                |
|  | LOT Batch code                      |
|  | Not made with natural rubber latex  |

**Rx ONLY**

Sterilized using E-Beam

Caution: Federal law (USA) restricts this device to sale or on the order of a physician.

**PRODUCT DESCRIPTION**

HaloGUARD™ Protective Disc with CHG is a sterile, single use disposable disc infused with the antibacterial agent chlorhexidine gluconate (CHG).

**INDICATIONS FOR USE:**

HaloGUARD™ Protective Disc with CHG is intended to cover insertion sites on adult patients. Common applications include IV catheters, central venous lines, epidural catheters, PICCs, hemodialysis catheters, orthopedic pins, other intravascular catheters and percutaneous devices.

**In-vitro Performance Testing**

HaloGUARD™ Protective Disc with CHG prevents contamination of the dressing only by inhibiting bacterial growth within the dressing. In-vitro studies have shown that HaloGUARD™ Protective Disc with CHG is effective at inhibiting hours of exposure to the following pathogens:

- *Staphylococcus aureus (MRSA)*
- *Pseudomonas aeruginosa*
- *Staphylococcus epidermidis (MRSE)*
- *Candida albicans*
- *Candida tropicalis*
- *Vancomycin-Resistant Enterococci (VRE)*
- *Enterobacter cloacae*
- *Klebsiella pneumoniae*



## PROTECTIVE DISC WITH CHG

STERILE R


 ATTWILL Medical Solutions  
 925 Development Drive, Lodi, WI  
 53555 www.ATTWILLmedical.com

To order products (customer service) in USA, call 1-608-592-6925.  
 For technical or product support, or any concerns, call 1-608-592-6925.

©2020 ATTWILL Medical Solutions. All rights reserved. HaloGUARD is a trademark of ATTWILL Medical Solutions

PK1354 Rev 5

#### CONTRAINDICATIONS, WARNINGS AND PRECAUTIONS

- HaloGUARD™ Protective Disc with CHG is intended for prescription use by clinicians and healthcare practitioners in a healthcare setting, as a sterile, single use disposable device.
- Do not use HaloGUARD™ Protective Disc with CHG on premature infants. Use of similar products on premature infants has resulted in hypersensitivity reactions and necrosis of the skin.
- The safety and effectiveness of HaloGUARD™ Protective Disc with CHG has not been established in children under 22 years of age. Use with pediatric populations is not recommended.
- Although it is not likely that other drugs or substances will interact with CHG as found in HaloGUARD™ Protective Disc with CHG, the clinician should verify if any medicines, vitamins, herbal products or other substances used by the patient may interact with CHG before using HaloGUARD™ Protective Disc with CHG.
- For external use only. Do not allow this product to contact the eyes, ears, mouth, or mucous membranes.
- Do not use HaloGUARD™ Protective Disc with CHG directly over a burn injury or on patients with a known sensitivity to CHG. Hypersensitivity reactions associated with the topical use of CHG have been reported in several countries. The most serious reactions have included anaphylaxis.
- Do not use HaloGUARD™ Protective Disc with CHG if patient is allergic to CHG. CHG may cause severe allergic reactions including wheezing, difficulty breathing, hives, shock, rash and facial swelling.
- Patient should stop use of HaloGUARD™ Protective Disc with CHG if irritation, sensitization, or allergic reaction occurs as these may be signs of a serious condition. Adverse reactions to CHG such as dermatitis, hypersensitivity, and generalized allergic reactions are rare, but if any such reactions occur, discontinue use of the dressing immediately.
- HaloGUARD™ Protective Disc with CHG is not intended to treat infections. Clinical studies on CHG to evaluate reduction in infections have not been performed.
- HaloGUARD™ Protective Disc with CHG is not intended to be used as a treatment of percutaneous device-related infections.
- Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician.

- HaloGUARD™ Protective Disc with CHG is intended for single use, DO NOT REUSE as it may create risk of patient infection, loss of essential design and material attributes of the device, compromise structural integrity, and may lead to device failure, patient injury, illness or death.

#### STORAGE AND HANDLING INFORMATION:

- Do not re-sterilize HaloGUARD™ Protective Disc with CHG using radiation.
- HaloGUARD™ Protective Disc with CHG is supplied sterile. Each unit package contains a single disc for single use only.
- Do not use if the package is opened or damaged. Do not use if seal is broken or compromised.
- Store between 15°C and 30°C (59°F and 86°F) in its original packaging.
- After use, handle and dispose of all unused product and packaging in accordance with accepted medical practice and applicable local, state, and federal laws and regulations. Discard any open and unused portions of the device.
- NOTE: HaloGUARD™ Protective Disc with CHG may appear yellow in color. This coloration does not reduce the antibacterial efficacy of the dressing.

#### SUMMARY OF SAFETY AND EFFECTIVENESS

HaloGUARD™ Protective Disc with CHG has undergone applicable in vitro, biocompatibility, sterility, and performance testing for safety and effectiveness.

| AVERAGE CONCENTRATION PER DISC |                           |                 |
|--------------------------------|---------------------------|-----------------|
| Model Number                   | Description               | CHG Amount (mg) |
| HG0141                         | 1 in disc, 4 mm hole      | 90              |
| HG0171                         | 1 in disc, 7 mm hole      | 85              |
| HG0075151                      | 0.75 in disc, 1.5 mm hole | 51              |